



Installation Guide for Origin Lab v2025b


- Download the software from the CCIS page. → [\(Download link\)](#) ←
- Extract the .zip folder.
- Run 'setup.exe' as administrator.
- Select the first option as shown below.




- Enter Serial number as “**GF3S5-3089-7911897**” as shown below.




- Continue as default settings as shown below.



- Continue as shown below.



- Specify location to save the user data. Keep default settings if no change is required.



- After the installation has finished. Open the software and update the licence server IP as '10.2.60.16' as shown below.

